



BRAEL - Medical Equipment

ul. Moscickiego 1, PL-24110 Pulawy,
Poland, tel./fax 48-81-4731393
www.brael.com.pl brael@brael.com.pl



Fetal Doppler: BabyPulse

| | |
|-----------------------|---|
| Model | FD-10D |
| Intended use | Fetal heart beat detection |
| Operating mode | Pulse Doppler |
| Power source | 2x AA alkaline battery |
| Operation time | 16 hours |
| Audio output | 1 W |
| Controls | Volume control and ON/OFF |
| Probe | Designed for general purpose (since 9 weeks) |
| Ultrasound frequency | 2MHz |
| Display method | 3-digit by 7-segment |
| FHR measurement range | 50-210 bpm |
| Dimensions | approx. 153x66x32 mm |
| Weight | approx. 220 g |
| Accessories | Gel, carrying case , 2x AA alkaline battery |
| Features | Auto-volume control Auto and manual FHR count Removing contact noise Battery indicator Pulse lamp |
| Order information | FD-10D Doppler with alkaline batteries FD-10D-R Doppler with NiMH rechargeable batteries+recharger |

CE 0044

TÜV NORD

2007-09-17