



## **BRAEL - Medical Equipment**

ul. Moscickiego 1, PL-24110 Pulawy,  
Poland, tel./fax 48-81-4731393  
[www.brael.com.pl](http://www.brael.com.pl) [brael@brael.com.pl](mailto:brael@brael.com.pl)



### Fetal Doppler: BabyPulse

Model	FD-10D
Intended use	Fetal heart beat detection
Operating mode	Pulse Doppler
Power source	2x AA alkaline battery
Operation time	16 hours
Audio output	1 W
Controls	Volume control and ON/OFF
Probe	Designed for general purpose (since 9 weeks)
Ultrasound frequency	2MHz
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Dimensions	approx. 153x66x32 mm
Weight	approx. 220 g
Accessories	Gel, carrying case , 2x AA alkaline battery
Features	Auto-volume control Auto and manual FHR count Removing contact noise Battery indicator Pulse lamp
Order information	FD-10D Doppler with alkaline batteries FD-10D-R Doppler with NiMH rechargeable batteries+recharger