



BRAEL - Medical Equipment

ul. Moscickiego 1, PL-24110 Pulawy,
Poland, tel./fax 48-81-4731393
www.brael.com.pl brael@brael.com.pl



Fetal Doppler: Doptone

Model: UDT-20D

with probe BWP2MHz

Model	UDT-20D
Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	2x AA alkaline battery
Operation time	32 hours
Audio output	1 W
Controls	Volume control and ON/OFF
Probe	Waterproof IPX7 – designed for water assisted labor and general purpose (since 10 weeks)
Ultrasound frequency	2MHz
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm on request 50-240 bpm
Dimensions	approx. 153x66x32 mm
Weight	approx. 220 g
Accessories	Gel, carrying case , 2x AA alkaline battery
Features	Auto-volume control Volume memorize Auto and manual FHR count Removing contact noise Battery indicator Pulse lamp Volume Memorize function at unit activation Resistant to interference from mobile phones
Order information	UDT-20WD Doppler with alkaline batteries UDT-20WD-R Doppler with NiMH rechargeable batteries+recharger