



## **BRAEL - Medical Equipment**

ul. Moscickiego 1  
ph.: +48 81 565 70 77  
[www.brael.com.pl](http://www.brael.com.pl)

PL-24-110 Pulawy  
fax: +48 81 565 70 78  
[brael@brael.com.pl](mailto:brael@brael.com.pl)



### Fetal Doppler: UDT-300

Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	230V/50Hz AC or 230V/50Hz AC + built-in battery
Audio output	10W
Controls	ON/OFF with volume control
Outputs	3.5 mm stereo jack to connect the headset
Interchangeable probe	B9P2MHz – wide-angle 9-crystal probe optimized for the quickest detection (since 14 weeks)
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Dimensions	255x110x210 mm
Weight	1,5 kg (UDT-300 SD) 2,2 kg (UDT-300 SD-R)
Accessories	One probe, gel
Accessories on demand	Additional probe and headset
Features	ARC Feature – the best signal control Signal Quality Indicator Objective fetal heart rate calculations Resistant to interference from mobile phones Built-in battery, above 6 hours working time (UDT-300SD-R)
Product code:	UDT-300 SD (230V/50Hz AC) UDT-300 SD-R (230V/50Hz AC + built-in battery)

2016-06-15

