



BRAEL - Medical Equipment

ul. Moscickiego 1, PL-24110 Pulawy,
Poland, tel./fax 48-81-4731393
www.brael.com.pl brael@brael.com.pl



Fetal Doppler: UDT-300

Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	230V/50Hz AC
Audio output	10W
Controls	ON/OFF with volume control
Outputs	3.5 mm stereo jack to connect the headset
Interchangeable probe	B7P2MHz – wide-angle 7-crystal probe optimized for the quickest detection (since 10 weeks)
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Dimensions	260x200x100 mm
Weight	1.5 kg
Accessories	One probe, gel
Accessories on demand	Additional probe and headset
Features	ARC Feature – the best signal control Signal Quality Indicator Objective fetal heart rate calculations Resistant to interference from mobile phones
Product code:	UDT-300SD

2011-09-23