



BRAEL - Medical Equipment

ul. Moscickiego 1
ph.: +48 81 565 70 77
www.brael.com.pl

PL-24-110 Puławy
fax: +48 81 565 70 78
brael@brael.com.pl



Fetal Doppler : UDT-300 with record

Model	UDT-300 with Record
Intended use	Fetal heart beat detection, FHR and fetal movement record
Operating mode	Continuous Doppler
Power source	230V/50Hz AC
Audio output	10W
Controls	ON/OFF with volume control and Record ON/OFF
Outputs	3.5 mm stereo jack to connect the headset, printer socket, 3.5 mm mono jack to connect the fetal movement marker
Interchangeable probes:	B9P2MHz – wide-angle 9-crystal probe
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Printer	Roll paper: 110mm 1cm/min
Dimensions of main unit	255x110x210 mm
Weight	Approx. 1.5 kg
Accessories	B9P2MHz probe , printer Kyoline , roll of paper, fetal movement marker , belt, gel
Accessories on demand	Additional narrow probe BP2MHz – optimized for deep penetration (since 10 weeks) and headset
Features	ARC Feature – the best signal control Objective fetal heart rate calculations FHR and Fetal Movement record
Order information	UDT-300SD-REC Doppler with 9-crystal probe and with record