



BRAEL - Medical Equipment

ul. Moscickiego 1
ph.: +48 81 565 70 77
www.brael.com.pl

PL-24-110 Puławy
fax: +48 81 565 70 78
brael@brael.com.pl



Fetal Doppler: UDT-300

Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	230V/50Hz AC or 230V/50Hz AC + built-in battery
Audio output	10W
Controls	ON/OFF with volume control
Outputs	3.5 mm stereo jack to connect the headset
Interchangeable probe	BP2MHz – optimized for deep penetration (since 10 weeks)
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Dimensions (WxHxD)	255x110x210 mm
Weight	1,5 kg (UDT-300D) 2,2 kg (UDT-300D-R)
Accessories	One probe, gel
Accessories on demand	Headset
Features	ARC Feature – the best signal control Signal Quality Indicator Objective fetal heart rate calculations Resistant to interference from mobile phones Built-in battery, above 6 hours working time (UDT-300D-R)
Product code:	UDT-300 D (230V/50Hz AC) UDT-300 D-R (230V/50Hz AC + built-in battery)

2016-06-15

