



## **BRAEL - Medical Equipment**

ul. Moscickiego 1, PL-24110 Pulawy,  
Poland, tel./fax 48-81-4731393  
[www.brael.com.pl](http://www.brael.com.pl) [brael@brael.com.pl](mailto:brael@brael.com.pl)



### Fetal Doppler UDT-300 with probe BP2MHz

Model	UDT-300
Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	230V/50Hz AC
Audio output	10W
Controls	ON/OFF with volume control
Outputs	3.5 mm stereo jack to connect the headset
Interchangeable probe	BP2MHz – optimized for deep penetration (since 9 weeks)
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Dimensions	260x200x100 mm
Weight	1.5 kg
Accessories	One probe, gel
Accessories on demand	Headset
Features	ARC Feature – the best signal control Signal Quality Indicator Objective fetal heart rate calculations
Order information	UD-300D

2009-06-12