

BRAEL-Aparaty Medyczne
 PL-24110 Pulawy, ul. Moscickiego 1
 tel./fax 48-81-8876948
www.brael.com.pl
brael@brael.com.pl



Fetal Doppler : UDT-300

Model	UDT-300
Intended use	Fetal heart beat detection
Operating mode	Continuous Doppler
Power source	230V/50Hz AC
Audio output	10W
Controls	ON/OFF with volume control and Record ON/OFF
Outputs	3.5 mm stereo jack to connect the headset
Interchangeable probes:	BP2MHz – optimized for deep penetration (since 10 weeks) B7P2MHz – wide-angle 7-crystal probe optimized for the quickest detection (since 10 weeks)
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm
Ultrasound frequency	2MHz
Dimensions	260x200x100 mm
Weight	1.5 kg
Accessories	One probe, gel
Accessories on demand	Additional probe and headset
Features	ARC Feature – the best signal control Objective fetal heart rate calculations
Order information	UD-300D Doppler with narrow angle probe UD-300SD Doppler with 7-crystal probe

CE 0044

4.03.2007