



## **BRAEL - Medical Equipment**

ul. Moscickiego 1  
ph.: +48 81 565 70 77  
[www.brael.com.pl](http://www.brael.com.pl)

PL-24-110 Pulawy  
fax: +48 81 565 70 78  
[brael@brael.com.pl](mailto:brael@brael.com.pl)



### Fetal Doppler: BabyPulse Model: **FD-10D**

Model	FD-10D
Intended use	Fetal heart beat detection
Operating mode	Pulse Doppler
Power source	2x AA alkaline battery
Operation time	16 hours
Audio output	1 W
Controls	Volume control and ON/OFF
Probe	Designed for general purpose (since 10 weeks)
Ultrasound frequency	2MHz
Display method	3-digit by 7-segment
FHR measurement range	50-210 bpm 50-240 bpm (option)
Dimensions (WxHxL)	approx. 65x32x145 mm
Weight	approx. 250 g (including batteries)
Accessories	Gel, carrying case , 2x AA alkaline battery
Features	Auto-volume control Auto and manual FHR count Removing contact noise Battery indicator Pulse lamp Auto shut off: 10 min in use/3 min no signal Resistant to interference from mobile phones
Order information	FD-10D Doppler with alkaline batteries FD-10D-R Doppler with NiMH rechargeable batteries+recharger



2016-06-15